

Dräger BabyFlow® plus Consumables and Accessories

A non-invasive respiratory support system designed to maximize performance and patient comfort.



PERFORMANCE

The BabyFlow plus respiratory support system was designed based on the results of state-of-the-art technology (Computational Flow Analysis) to minimize the imposed work of breathing for the patient. The streamlined adapter head optimizes the flow distribution.

ENHANCED COMFORT FOR PATIENT AND CAREGIVER

We want your patients to feel comfortable. The less discomfort a patient experiences, the greater the opportunity for healthy growth and development during this critical time in the patient's recovery.

The adjustable angle allows you to fit the prongs and masks properly on the patient, thus reducing leaks and providing more stable CPAP levels. More stable CPAP levels reduce noise levels, making BabyFlow plus a low noise level respiratory support system, therefore, supporting Developmental Care practices and offering the best possible therapy for tiny patients.*

EASE OF USE

The ease of use of BabyFlow plus in combination with innovative components offer effective and gentle therapy where it's needed.

The BabyFlow plus respiratory support system is designed to address the current issues with non-invasive respiratory support therapy. Every single component, including individually sized masks and nasal prongs, has been engineered to maximize performance and patient comfort.

The flexible, kink resistant hoses enable optimal positioning of the patient. Swivel connectors reduce the necessary disconnections of the breathing circuit to a minimum, thus ensuring stable CPAP pressures.

In combination with the Dräger Babylog® VN500, the BabyFlow plus supports a comprehensive selection of modalities for non-invasive respiratory support. Babylog VN500 automatically adjusts the flow according to the existing leaks to ensure stable CPAP pressure and prevent nuisance alarms.

- * L. Kirchner, M. Wald, V. Jeitler, A. Pollack, "In vitro Comparison of Noise Levels produced by different CPAP generators", in: Neonatology, 2012; 101: 95-100
- * Karam et al., "Noise levels during nCPAP are flow-dependent but not device dependent", in: Archives of Disease in Childhood Neonatal Edition, 2008; 93: F 132-F134



BabyFlow® for disposable use with measuring tape

ORDERING INFORMATION

BabyFlow® plus disposable

S	20 pieces	MP03700
M	20 pieces	MP03701
L	20 pieces	MP03702

Microfiber Caps

XS	white	5 pieces	84 18 533
S	yellow seam	5 pieces	84 18 534
M	red seam	5 pieces	84 18 535
L	green seam	5 pieces	84 18 536
XL	orange seam	5 pieces	84 18 537
XXL	blue seam	5 pieces	84 18 538
XXL	plus white/black seam	5 pieces	84 18 539

Fixation Straps

Fixation Strap	13 cm	10 pieces	MP03870
Fixation Strap	15 cm	10 pieces	MP03871
Fixation Strap	18 cm	10 pieces	MP03872
Fixation Strap	20 cm	10 pieces	MP03873

Headgear

XXS	5 pieces	MP03880
XS	5 pieces	MP03881
S	5 pieces	MP03882
M	5 pieces	MP03883
L	5 pieces	MP03884
XL	5 pieces	MP03885

Neo mask

S	10 pieces	84 18 491
M	10 pieces	84 18 490
L	10 pieces	84 18 619
4	10 pieces	MP01494
5	10 pieces	MP01495
6	10 pieces	MP01496

Standard prongs

Prongs XS	10 pieces	84 18 415
Prongs S	10 pieces	84 18 605
Prongs M	10 pieces	84 18 416
Prongs L	10 pieces	84 18 531
Prongs XL	10 pieces	84 18 417
Prongs XXL	10 pieces	84 18 617

Wide prongs

Prongs 3	10 pieces	84 18 603
Prongs 4	10 pieces	84 18 607
Prongs 5	10 pieces	84 18 609
Prongs 6	10 pieces	84 18 611
Soft connector N	20 pieces	MP03826

D-7478-2016

**BabyFlow® plus disposable**

D-99280-2013

**Microfiber caps**


D-23988-2015

**Headgear**

D-23992-2015/D34333-2009

**Neo masks**

MT-2010-2005

**Prongs**

Demo pack BabyFlow® plus, size S**MP03703**

BabyFlow® plus, size S, 10 pieces

2 x Prongs XS

2 x Prongs S

2 x Prongs M

2 x Prongs 3

2 x Neo mask S

2 x Neo mask M

2 x Neo mask L

2 x Headgear XXS

2 x Headgear XS

2 x Headgear S

Tape measure and sizing gauge

Soft connector N, 20 pieces

Demo pack BabyFlow® plus, size M**MP03704**

BabyFlow® plus, size M, 10 pieces

2 x Prongs M

2 x Prongs L

2 x Prongs XL

2 x Prongs 4

2 x Neo mask M

2 x Neo mask L

2 x Neo mask 4

2 x Headgear S

2 x Headgear M

2 x Headgear L

Tape measure and sizing gauge

Soft connector N, 20 pieces

Demo pack BabyFlow® plus, size L**MP03705**

BabyFlow® plus, size L, 10 pieces

2 x Prongs XL

2 x Prongs XXL

2 x Prongs 5

2 x Prongs 6

2 x Neo mask 4

2 x Neo mask 5

2 x Neo mask 6


2 x Headgear M

2 x Headgear L

2 x Headgear XL

Tape measure and sizing gauge

Soft connector N, 20 pieces

**BabyFlow® plus Demo Pack****BabyFlow® plus complete solution**

D-48840-2015



Application of BabyFlow plus with mask on full-term infant

D-48856-2015



Application of BabyFlow plus on pre-term infant

Not all products, features, or services are for sale in all countries.
Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to www.draeger.com/trademarks to find the current status.

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

USA

Draeger, Inc.
3135 Quarry Road
Telford, PA 18969-1042
Tel +1 800 4DRAGER
(+1 800 437 2437)
Fax +1 215 723 5935
info.usa@draeger.com

Manufacturer:

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany