




Conox®

Smart Anaesthesia Monitoring

Data sheet

Two monitoring solutions in one device

Advanced digital processing algorithms
for depth of anaesthesia monitoring

Both hypnotic effect and probability of
response indices are measured from the
same electrodes



Monitor main characteristics

Dimensions

215 x 150 x 65 mm

Weight

785 g

Display size

7" (154.08 x 85.92 mm)

Battery

3.7 Volts Li-ion battery
Capacity 3500 mAh
Rechargeable

Bluetooth®

2.1, selectable, send only, up to 8 m

Patient Cable length

2 m

Technical specifications - EEG acquisition

EEG resolution

44.58 nV

Input noise

≤ 2.9 µV peak to valley

CMRR

89 dB

Sampling frequency

1024 Hz

Input signal range

+/- 374 mV

Resolution

24 bit ADC

Indices performance

qCON index range

0 - 99

qNOX index range

0 - 99

EMG index range

0 - 100

SQL index range

0 - 100

BSR index range

0 - 100

Indices update rate

Every 1 s

qCON index alarms

Selectable, high and low limit, audio and visual

qCON / qNOX indices update time

Current values correspond to EEG signal down to 20 s, i.e. displayed indices represent a moving average of 20 s of patient data.

Indices trend graphing

Selectable, qCON, qNOX, EMG, BSR

Lead off detection

Continuous

Sensor impedance check

Automatic every 15 min or on request via panel button

Conox and Conox logo are registered trademarks of Fresenius Kabi. Due to our policy of continuous product and specification development as well as changes in the regulation standards, the features described above are subject to change. Please consult our Sales Executive for the most updated information.



**FRESENIUS
KABI**
caring for life



Fresenius Kabi AG
61346 Bad Homburg, Germany
Tel: +49 (0) 61 72 / 686-0
www.fresenius-kabi.com



Quantum Medical S.L.U.
Av Ernest Lluch 32, Tower 2
Tecnocampus, 3rd Floor, Office 317,
08302 Mataró, Spain
Tel: +34 93 702 1950
www.quantummedical.com

